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
Dr. Thomas Härtling
Telefon 0351 88815-550
thomas.haertling@ikts.fraunhofer.de

www.ikts.fraunhofer.de

FRAUNHOFER INSTITUTE FOR CERAMIC TECHNOLOGIES AND SYSTEMS IKTS

1 Medical product wrapped with a foil containing ceramic pigments. The green luminescence of the latter is clearly visible.

2, 3 X-Ray microscopic pictures showing the distribution of ceramic pigments inside the packaging material (a polymer foil).


STERILIZATION PROCESS MONITORING

Sterilization based on the impact of low-energy electrons is a successful alternative to current techniques for sterilizing highly sensitive medical products, thermo sensitive polymers, electronics, or functional biological materials. Fraunhofer IKTS is now presenting a fast and secure procedure for monitoring the electron beam (e-beam) assisted sterilization of such objects.

Approach

Core elements of the monitoring system are ceramic pigments integrated into or onto the packaging material, e.g. a polymer foil. These pigments change their optical properties upon e-beam treatment which allows a derivation of the applied electron dose. The latter can be read out either by means of a hand-held device or an integrated optical system (Figure 1).

Working principle based on the change of the luminescence decay time of ceramic pigments

Advantages

In contrast to standard techniques like film dosimeters, no time-consuming analysis is required and results are available immediately after the sterilization process. Hence, an inline integration of the control system into the e-beam plant is feasible. Even a position-dependent analysis of e-beam sterilization is possible. Thus, a local sterilization of previously omitted surface areas can be realized.

Fraunhofer Institute for Ceramic Technologies and Systems IKTS

Maria-Reiche-Strasse 2
01109 Dresden, Germany

Contact

Dr. Thomas Härtling
Phone +49 351 88815-550
thomas.haertling@ikts.fraunhofer.de

www.ikts.fraunhofer.de